

## Xpeedy™ DOA 20

<b>Reagent Dispensing</b>	Pre-mixing and dispensing	
<b>Robotic Arms</b>	2-Channel Capper/Decapper and 2-Tip Dispenser	
	Pressure Liquid Surface Detection	
	Dispensing range: 5-1000 uL	
	Dispensing Coefficient of Variance (CV)	Dispensing Accuracy (A)
	20uL: ≤2%	20uL: ±3%
	100uL: ≤1.25%	100uL: ±2%
	500uL: ≤0.4%	500uL: ±1%
	1000uL: ≤0.4%	1000uL: ±1%
<b>Mixing Module</b>	3 speed options	
<b>Heating Module</b>	Temperature Control	Room temperature +5 °C ~ 105 °C
	Temperature Accuracy	≤±0.5 °C
<b>Output</b>	LCD Touch-screen	
<b>Connectivity</b>	USB Port	
<b>Power</b>	100-240 V AC, 50/60 Hz, 400 W	
<b>Dimension/Weight</b>	112.6 cm x 77.5 cm x 80 cm/176 kg	