


# REVOLUTIONIZING HEALTHCARE

ONE AUTOMATED TESTING LAB AT A TIME





## ABOUT US

---

CURA Diagnostics is a groundbreaking in-vitro diagnostics company committed to pushing the boundaries of healthcare efficiency and accessibility. Through the development and provision of a diversified portfolio of first-in-class or best-in-class diagnostic tests, our company is at the forefront of bringing innovative and transformative solutions to the medical diagnostic sector. With core commitments to be **C**ustomer-centered, **U**ser-friendly, **R**eliable, and **A**ffordable, CURA Diagnostics aims to transform the standards of care through the power of automation and cutting-edge technologies.

CURA's automated lab testing systems will vastly increase efficiency and productivity, reduce turn-around time, as well as reduce operation costs enabling quick, accurate diagnoses for effective treatment planning. CURA Diagnostics further provides a range of turn-key service solutions for laboratory startup, validation, and compliance management.

Through our products and services, we aim to revolutionize the way healthcare providers identify, diagnose, and treat illnesses, improving the quality of life for millions of patients worldwide.





## Xpeedy™ DOA 20

Automated Sample Preparation for  
Triple Quadrupole LC/MS/MS Analysis

### Technical Specification

• Reagent Dispensing	Pre-mixing and Dispensing	
• Robotic Arms	2-Channel Capper/Decapper and 2-Tip Dispenser	
	Pressure Liquid Surface Detection	
	Dispensing Range: 5-1000 µL	
	Dispensing Coefficient of Variance (CV)	Dispensing Accuracy (A)
• Liquid Dispenser	20 µL: ≤ 2%	20 µL: ± 3%
	100 µL: ≤ 1.25%	100 µL: ± 2%
	500 µL: ≤ 0.4%	500 µL: ± 1%
	1000 µL: ≤ 0.4%	1000 µL: ± 1%
• Mixing Module	3 Speed Options	
	Temperature Control	Room temperature +5 °C ~ 105 °C
• Heating Module	Temperature Accuracy	≤ ± 0.5 °C
• Output	LCD Touch-screen	
• Connectivity	USB Port	
• Power	100-240 V AC, 50/60 Hz, 400 W	
• Dimension/Weight	112.6 cm x 77.5 cm x 80 cm/176 kg	

# CLINICAL LAB SERVICES



## Start

Laboratory Startup, CLIA Licensure, Equipment, Consumables Supplies, and Lab Design Solutions.



## Run

Lab and Quality Management, Accreditation Requirements, Test Validation, and Business Intelligence Services.



## Grow

Customization, Adding Lines of Testing, Laboratory Management, and Commercialization Solutions.



+1-339-999-2959



[www.cura-dx.com](http://www.cura-dx.com)



[info@cura-dx.com](mailto:info@cura-dx.com)



800 W Cummings Park,  
Suite 3550, Woburn,  
MA 01801